




# STARCUT

Semi-automatic biopsy needle STARCUT



Quality, All the way, Always, That's TSK.



# STARCUT

## Semi-automatic biopsy needle STARCUT

The semi-automatic biopsy needle STARCUT is used to insert into the human body to collect soft tissue for the purpose of examination, treatment, or diagnosis.

### Feature

#### 1 Needle tip

Sharp needle tip reduces stress during puncture.

#### 2 Depth Marker


You can visually check the puncture depth. Lines are placed at 1cm intervals, and double lines are placed every 5cm.

#### 3 Two different penetration lengths

The penetration length can be selected from 11 mm and 22 mm to suit the procedure.

#### 4 Lock position

The set protrusion length of 11 mm and 22mm can be visually confirmed from the window of the main unit. (写真は22mm)



※Please check the latest package inserts before use.

### Product standard

Size		Product code
		STARCUT
14G	90mm	SAG-14090
	130mm	SAG-14130
	160mm	SAG-14160
	220mm	SAG-14220
16G	90mm	SAG-16090
	130mm	SAG-16130
	160mm	SAG-16160
	220mm	SAG-16220
18G	90mm	SAG-18090
	130mm	SAG-18130
	160mm	SAG-18160
	220mm	SAG-18220
20G	90mm	SAG-20090
	130mm	SAG-20130
	160mm	SAG-20160
	220mm	SAG-20220

#### TSK Guide Needle for Semi-Automatic Biopsy Needle STARCUT

Echo-compatible coaxial guide needle.

Please inquire separately for size and price.



Gamma sterilization 10 pcs/box



Manufactured and Distributed by TSK

1510-1 Soja-machi, Tochigi-shi, Tochigi-Ken, 328-0002 Japan

TEL: +81 282 27 8426 FAX: +81 282 25 6511

Distributed by :